

# PCR FOR GASTROINTESTINAL VIRAL PATHOGEN PANEL

(61028162, PCRG12)

## SPECIMEN & CONTAINER

### Stool

Stool container  
(without preservative)

## TRANSPORTATION REQUIREMENT

2 – 8°C (up to 2 days)  
– 20°C (up to 1 month)

## METHOD

Multiplex Real-Time PCR

## TEST SCHEDULE

Daily (Monday to Saturday)  
Cut off time: 11.00 AM

## TURNAROUND TIME (TAT)

1 – 2 Days  
upon specimen arrival

## DETECTION OF

1. Norovirus G1 (NVG1)
2. Norovirus G2 (NVG2)
3. Rotavirus A (ROV)
4. Adenovirus 40/41 (ADV-F)
5. Astrovirus (ASV)
6. Sapovirus (SV)