



## ENDOCRINE AUTOIMMUNE PROFILE

### DIABETES MELLITUS AUTOIMMUNE PROFILE 1 (AiDMp1, 61050917)

#### TEST LIST

- Islet cell cytoplasmic Autoantibody, ICA
- Glutamic acid decarboxylase Autoantibody, GADA
- Insulinoma associated protein-2 Autoantibody, IA-2A
- Insulin Autoantibody, IAA

#### METHOD

ELISA  
ELISA  
ELISA  
ELISA

#### SPECIMEN

2 mL of separated serum. Refrigerated (2-8°C) or freeze (-20°C) immediately.

#### TRANSPORTATION REQUIREMENT

Refrigerated (for Klang Valley customer) or frozen (for outstation customer).

Sample shall arrive Lablink Central within 24 hours.

#### TURNAROUND TIME (TAT)

2 – 3 working days upon receiving sample at Lablink Central.



## ENDOCRINE AUTOIMMUNE PROFILE

### THYROID AUTOANTIBODY SCREEN PROFILE (AiThyP1, 61050941)

TEST LIST	METHOD
• Thyrotropin Receptor Antibody (TSH Receptor Antibody - TSHRepAb)	CLIA
• Thyroid Peroxidase (Anti-Thyroid peroxidase - Anti-TPO/AMC)	CMIA
• Thyroglobulin (Anti-Thyroglobulin - Anti-TG)	CMIA

#### SPECIMEN

2 mL of separated serum. Refrigerated (2-8°C) or freeze (-20°C) immediately.

#### TRANSPORTATION REQUIREMENT

Refrigerated (for Klang Valley customer) or frozen (for outstation customer).

Sample shall arrive Lablink Central within 24 hours.

#### TURNAROUND TIME (TAT)

2 – 3 working days upon receiving sample at Lablink Central.