



# PHADIATOP ALLERGEN SCREEN

(S, 61040301)

## ALLERGEN NAME

Phadiatop Allergen

Qualitative screening tests useful for identifying patients with a high likelihood of an allergic disease (so called atopic individuals).

Phadiatop utilize balanced mixtures of relevant allergens, often implicated in causing allergic diseases in adults and in young children, respectively.

*Ref: Reference: Williams PB, et al. Ann Allergy Asthma Immunol. 2001;86:196-202*

## SPECIMEN

Minimum 1ml Serum - Refrigerated (2-8°C) or freeze (-20°C) immediately

## TEST SCHEDULE

By Batch: Wednesday (Before 11 a.m)  
Report on Thursday

## TURNAROUND TIME (TAT)

7 working days

## TRANSPORTATION REQUIREMENT

Ship refrigerated (for Klang Valley) or frozen (compulsory for outstation).  
Sample to arrive Lablink Central within 24 hours.

## METHOD

ImmunoCAP EIA



"1 out of 3 Malaysians is currently suffering from some form of allergy".  
Assoc Prof Dr Ranbir Kaulsay (President of Malaysian Society of Allergy & Immunology, MSAI)



# DOMESTIC ALLERGY PANEL

(GP94D, 61050636)

## ALLERGEN NAME

### TOTAL IgE

### INHALATION

<b>Mites</b>	Dermatophagoides pter. Dermatophagoides farinae. House dust Blomia tropicalis
<b>Animals/ Insects</b>	Cockroach

## SPECIMEN

Minimum 2ml Serum - Refrigerated (2-8°C) or freeze (-20°C) immediately

## TEST SCHEDULE

By Batch: Wednesday (Before 11 a.m)  
Report on Thursday

## TURNAROUND TIME (TAT)

7 working days

## TRANSPORTATION REQUIREMENT

Ship refrigerated (for Klang Valley) or frozen (compulsory for outstation).  
Sample to arrive Lablink Central within 24 hours.

## METHOD

ImmunoCAP EIA

"1 out of 3 Malaysians is currently suffering from some form of allergy " .  
*Assoc Prof Dr Ranbir Kaulsay (President of Malaysian Society of Allergy & Immunology, MSAI)*



# SEAFOOD FOOD ALLERGY PANEL

(GP94E, 61050637)

## ALLERGEN NAME

### TOTAL IgE

### SEA FOOD

<b>Fish</b>	Anchovy
<b>Shellfish</b>	Crab Shrimp Clam
<b>Squid</b>	

## SPECIMEN

Minimum 2ml Serum - Refrigerated (2-8°C) or freeze (-20°C) immediately

## TEST SCHEDULE

By Batch: Wednesday (Before 11 a.m)  
Report on Thursday

## TURNAROUND TIME (TAT)

7 working days

## TRANSPORTATION REQUIREMENT

Ship refrigerated (for Klang Valley) or frozen (compulsory for outstation).  
Sample to arrive Lablink Central within 24 hours.

## METHOD

ImmunoCAP EIA

"1 out of 3 Malaysians is currently suffering from some form of allergy " .  
*Assoc Prof Dr Ranbir Kaulsay (President of Malaysian Society of Allergy & Immunology, MSAI)*



# PEDIATRIC FOOD ALLERGY PANEL

(GP94F, 61050638)

## ALLERGEN NAME

### TOTAL IgE

### SEA FOOD

**Fish**

Anchovy

### FOOD

**Nuts**

Peanut  
Soybean

**Meats /  
Poultry**

Chicken

**Egg/Milk  
Yeast**

Egg white  
Egg yolk  
Cow's Milk

**Protein/  
Grains/  
Additiv./  
Cooking Oil**

Wheat

## SPECIMEN

Minimum 2ml Serum - Refrigerated (2-8°C) or freeze (-20°C) immediately

## TEST SCHEDULE

By Batch: Wednesday (Before 11 a.m)  
Report on Thursday

## TURNAROUND TIME (TAT)

7 working days

## TRANSPORTATION REQUIREMENT

Ship refrigerated (for Klang Valley) or frozen (compulsory for outstation).  
Sample to arrive Lablink Central within 24 hours.

## METHOD

ImmunoCAP EIA

"1 out of 3 Malaysians is currently suffering from some form of allergy ".  
Assoc Prof Dr Ranbir Kaulsay (President of Malaysian Society of Allergy & Immunology, MSAI)



# MALAYSIAN ALLERGY PANEL

(GP94G, 61050651)

## ALLERGEN NAME

INHALATION		FOOD		SEA FOOD	
<b>Grasses</b>	Bermuda grass Johnson grass	<b>Nuts</b>	Peanut Soybean	<b>Fish</b>	Tuna Anchovy
<b>Tree</b>	Rubber Tree (Latex)	<b>Meats / Poultry</b>	Beef Chicken		
<b>Mites</b>	Dermatophagoides pter. Dermatophagoides farinae. House dust	<b>Egg/Milk Yeast</b>	Egg white Egg yolk Cow's Milk		
<b>Animals/ Insects</b>	Cat Dog Cockroach	<b>Protein/ Grains/ Additiv./ Cooking Oil</b>	Wheat	<b>Shellfish</b>	Crab Shrimp Clam
<b>Molds</b>	Candida albicans Cladosporium herbarum Aspergillus fumigatus Alternaria alternate Mucor Racemosus Penicillium B	<b>Vege. / Fruits</b>	Banana Orange		

## SPECIMEN

Minimum 2ml Serum - Refrigerated (2-8°C) or freeze (-20°C) immediately

## TEST SCHEDULE

By Batch: Wednesday (Before 11 a.m)  
Report on Thursday

## TURNAROUND TIME (TAT)

7 working days

## TRANSPORTATION REQUIREMENT

Ship refrigerated (for Klang Valley) or frozen (compulsory for outstation).  
Sample to arrive Lablink Central within 24 hours.

## METHOD

ImmunoCAP EIA

"1 out of 3 Malaysians is currently suffering from some form of allergy ".  
Assoc Prof Dr Ranbir Kaulsay (President of Malaysian Society of Allergy & Immunology, MSAI)



# ATOPY ALLERGY PANEL

(GP94H,61050887)

## ALLERGEN NAME

INHALATION		FOOD		SEA FOOD	
<b>Grasses</b>	Bermuda grass Timothy grass Grass mix 5: Sweet vernal grass Bermuda grass Timothy grass Cultivated rye	<b>Nuts</b>	Peanut Soybean Hazelnut Almond	<b>Fish</b>	Codfish Tuna Salmon
<b>Tree</b>	Acacia Pine (Australia) Kapak	<b>Meats / Poultry</b>	Duck meat Beef (cooked) Pork (cooked) Chicken Lamb	<b>Shellfish</b>	Crab Shrimp (Pacific) Lobster Shellfish mix 1 : Spiny lobster Oyster Clam
<b>Mites</b>	Dermatophagoides pter. Dermatophagoides farinae. Dermatophagoides microceras Tyrophagus putrescentiae Glycyphagus domesticus Blomia tropicalis	<b>Egg/Milk Yeast</b>	Egg white Cow's Milk nBos d4 alpha-lactalbumin (Milk) nBos d5 beta-lactoglobulin (Milk) nBos d8 Casein (Milk) nBos d6 BSA (Milk) Cheddar Cheese Baker's yeast	<b>INSECT VENOMS</b>	
<b>Animals/ Insects</b>	Cat Dog Horse Feather mix 1: Chicken feathers, Duck feathers, Goose feathers Coackroach (German)	<b>Protein/ Grains/ Additiv./ Cooking Oil</b>	Wheat flour Gluten Chocolate Glutamate Coffee Oil palm	<b>Insect</b>	Honey bee
<b>Molds</b>	Candida albicans Mould mix 1: Penicillium notatum Cladosporium herbarum Aspergillus fumigatus Alternaria alternata	<b>Vege. / Fruits</b>	Tomato Garlic Strawberry Kiwi	<b>CROSS REACTIVITY</b>	
				<b>Marker</b>	CCD Marker

## SPECIMEN

Minimum 1ml Serum - Refrigerated (2-8°C) or freeze (-20°C) immediately

## TEST SCHEDULE

By Batch: Wednesday (Before 11 a.m)  
Report on Thursday

## TURNAROUND TIME (TAT)

7 working days

## TRANSPORTATION REQUIREMENT

Ship refrigerated (for Klang Valley) or frozen (compulsory for outstation).  
Sample to arrive Lablink Central within 24 hours.

## METHOD

Immunoblot