



VASCULITIS & IMMUNE NEPHRITIS AUTOIMMUNE PROFILES

VASCULITIS AUTOIMMUNE PROFILE 1 (AiVasP1, 61050923)

TEST LIST

- ANCAIF - Anti-neutrophil cytoplasmic antibody
- Autoimmune Vasculitis Screening:
 - ANCApMPO - Myeloperoxidase (MPO) pANCA antibody
 - ANCAcPR3 - Proteinase 3 (PR3) cANCA antibody
 - GBM - Glomerular basement membrane (GBM) antibody

METHOD

IFA/IFT

ELISA

ELISA

ELISA

SPECIMEN

2 ml of separated serum. Refrigerated (2-8°C) or freeze (-20°C) immediately.

TRANSPORTATION REQUIREMENT

Refrigerated (for Klang Valley customer) or frozen (for outstation customer).
Sample shall arrive Lablink Central within 24 hours.

TURNAROUND TIME (TAT)

2 – 3 working days upon receiving sample at Lablink Central.



VASCULITIS & IMMUNE NEPHRITIS AUTOIMMUNE PROFILES

VASCULITIS AUTOIMMUNE PROFILE 2 (AiVasP2, 61050924)

TEST LIST

- ANAIF - Anti-nuclear antibody (ANA) IFA/IIFT
- ANCAIF - Anti-neutrophil cytoplasmic antibody, ANCA IFA/IIFT
- Autoimmune Vasculitis Screening:
 - ANCApMPO - Myeloperoxidase (MPO) pANCA antibody
 - ANCAcPR3 - Proteinase 3 (PR3) cANCA antibody
 - GBM - Glomerular basement membrane (GBM) antibody

METHOD

IFA/IIFT
IFA/IIFT
ELISA
ELISA
ELISA

SPECIMEN

2 ml of separated serum. Refrigerated (2-8°C) or freeze (-20°C) immediately.

TRANSPORTATION REQUIREMENT

Refrigerated (for Klang Valley customer) or frozen (for outstation customer).
Sample shall arrive Lablink Central within 24 hours.

TURNAROUND TIME (TAT)

2 – 3 working days upon receiving sample at Lablink Central.

Additional 2 – 3 working days for Positive Anti-nuclear antibody (ANA) IFA/IIFT for titration procedure.



VASCULITIS & IMMUNE NEPHRITIS AUTOIMMUNE PROFILES

VASCULITIS AUTOIMMUNE PROFILE (AiVasP3, 61050925)

TEST LIST	METHOD
• ANAIF - Anti-nuclear antibody (ANA) IFA/IIFT	IFA/IIFT
• ANCAIF - Anti-neutrophil cytoplasmic antibody, ANCA IFA/IIFT	IFA/IIFT

SPECIMEN

1 ml of separated serum. Refrigerated (2-8°C) or freeze (-20°C) immediately.

TRANSPORTATION REQUIREMENT

Refrigerated (for Klang Valley customer) or frozen (for outstation customer).
Sample shall arrive Lablink Central within 24 hours.

TURNAROUND TIME (TAT)

2 – 3 working days upon receiving sample at Lablink Central.

Additional 2 – 3 working days for Positive Anti-nuclear antibody (ANA) IFA/IIFT for titration procedure.



VASCULITIS & IMMUNE NEPHRITIS AUTOIMMUNE PROFILES

VASCULITIS AUTOIMMUNE PROFILE 4 (AiVasP4, 61050926)

TEST LIST

- ANCApMPO - Myeloperoxidase (MPO) pANCA antibody
- ANCAcPR3 - Proteinase 3 (PR3) cANCA antibody
- GBM - Glomerular basement membrane (GBM) antibody

METHOD

ELISA
ELISA
ELISA

SPECIMEN

1 ml of separated serum. Refrigerated (2-8°C) or freeze (-20°C) immediately.

TRANSPORTATION REQUIREMENT

Refrigerated (for Klang Valley customer) or frozen (for outstation customer).
Sample shall arrive Lablink Central within 24 hours.

TURNAROUND TIME (TAT)

2 – 3 working days upon receiving sample at Lablink Central.



VASCULITIS & IMMUNE NEPHRITIS AUTOIMMUNE PROFILES

NEPHRITIS AUTOIMMUNE PROFILE 5 (AiVasP5, 61050927)

TEST LIST	METHOD
• ANAIF - Anti-nuclear Antibody (ANA) IFA/IIFT	IFA/IIFT
• ANCAIF - Anti-neutrophil cytoplasmic Antibody, ANCA IFA/IIFT	IFA/IIFT
• DNAAb - Ds-DNA Antibody	ELISA
• Autoimmune Vasculitis Screening:	
• ANCApMPO - Myeloperoxidase (MPO) pANCA Antibody	ELISA
• ANCAcPR3 - Proteinase 3 (PR3) cANCA Antibody	ELISA
• GBM - Glomerular basement membrane (GBM) Antibody	ELISA
• Autoimmune Renal Membranous Nephropathy Screening:	
• Anti-phospholipase-A2-receptor (PLA2R) Antibody	IFA/IIFT

SPECIMEN

1 ml of separated serum. Refrigerated (2-8°C) or freeze (-20°C) immediately.

TRANSPORTATION REQUIREMENT

Refrigerated (for Klang Valley customer) or frozen (for outstation customer).
Sample shall arrive Lablink Central within 24 hours.

TURNAROUND TIME (TAT)

2 – 3 working days upon receiving sample at Lablink Central.
Additional 2 – 3 working days for Positive Anti-nuclear antibody (ANA) IFA/IIFT for titration procedure.