

# PCR FOR GASTROINTESTINAL BACTERIAL PATHOGEN PANEL

(61028163, PCRG13)

## SPECIMEN & CONTAINER

### Stool

Stool container  
(without preservative)

## TRANSPORTATION

### REQUIREMENT

2 – 8°C (up to 2 days)  
– 20°C (up to 1 month)

## METHOD

Multiplex Real-Time PCR

## TEST SCHEDULE

Daily (Monday to Saturday)  
Cut off time: 11.00 AM

## TURNAROUND TIME (TAT)

1 – 2 Days  
upon specimen arrival

## DETECTION OF

1. *Campylobacter spp. (Cam)*
2. *Clostridium difficile toxin B (CdB)*
3. *Salmonella spp. (Sal)*
4. *Shigella spp./EISC (Sh/EI)*
5. *Vibrio spp. (Vib)*
6. *Yersinia enterocolitica (Yer)*
7. *Aeromonas spp. (Aer)*