



**MEMORANDUM**

To : All Valued Clients of Lablink Medical Laboratory  
 From : Pathologist In-Charge  
 Date : 15<sup>th</sup> January 2018  
 Reference : LABLINK/PRMCS/MEMO/2018-005

**SUBJECT : INTRODUCTION OF NEW AUTOIMMUNE TESTS**

The above subject is kindly referred.

With effective of 15<sup>th</sup> January 2018, Lablink Medical Laboratory is pleased to introduce new autoimmune tests. Details of the new autoimmune profiles are as follows:-

NO	PROFILE	TEST CODE		
1	<b>Connective Tissue Diseases &amp; Rheumatology Autoimmune Profiles</b>	- AiCTDp1	- AiCTDp4	- AiCTDp5
2	<b>Vasculitis &amp; Immune Nephritis Autoimmune Profiles</b>	- AiVasP1 - AiVasP2	- AiVasP3 - AiVasP4	- AiVasP5
3	<b>Neurology &amp; Encephalitis Autoimmune Profiles</b>	- AiNeuroP1 - AiNeuroP2 - AiNeuroP3	- AiNeuroP4 - AiNeuroP5 - AiNeuroP6	- AiNeuroP7 - AiNeuroP8 - AiNeuroP9
4	<b>Endocrine Autoimmune Profile</b>	-AiDMp1		
5	<b>Gastroenterology Autoimmune Profiles</b>	- AiGasP1	- AiGasP2	- AiGasP3
6	<b>Individual Test</b>	- ANAIF - DNAAb - ANCAIF - NMDARcsf - IA2A	- NMDARab - AQP4 - ANCApMPO - ANCAcPR3 - InsulinAb	- GBM - PLA2RAb - IsletAb - GADA - AChRAb

To ensure rapid diagnosis, Lablink promises a laboratory report ready within 72 hours once specimen received in Lablink Central. For a positive sample of IFA/IIFT, 72 hours additional Turn-Around Time (TAT) is required.

Kindly refer to the attachment for details on the method, specimen requirement and price of test. Any questions regarding the test should be directed to the on-site Laboratory or marketing services at 03-4023 4588.

Thank you.

MUHAMMAD NAZRI AZIZ